



Manufacture's Declaration of Laser Radiation of Safe Laser 580 Duo

Company: **Laser Light Medical Innovation Ltd.**

Address: **Hungary, 1138 Budapest, Hajóépítő sétány 4. B.IX.201.**

as manufacturer we declare that our product

Name: **Veterinary therapy laser device**
Safe Laser 580 Duo

is suitable for professional and for personal use.

It radiates special scattered laser light which is safe for human eyes according to the standard:
Safety of laser products - Part 1: Equipment classification and requirements (IEC 60825-1:2014)

Laser safety: Class 3/R

Some countries have import restrictions on hand held, battery operated, laser pointers with an assessable emission level greater than 1mW, which emits dangerous parallel laser beam.

The Safe Laser 580 Duo device is NOT a “laser pointer”.

Hungary, Budapest, 11 July 2025

A handwritten signature in blue ink, appearing to read "Rozsa Tamas", with a stylized flourish at the end.

Tamas Rozsa

CEO Laser Light Medical Innovation Ltd.

For more information contact: info@safelaser.eu